Birmingham Knee Replacement

EN

INSTRUCTIONS FOR USE

Carefully read all instructions and be familiar with the surgical technique prior to use.

Description and Intended Use

The Birmingham Knee Replacement (BKR) is a rotating platform, posterior stabilised, primary total knee prosthesis, for use when both cruciate ligaments are excised.

The BKR consists of femoral, tibial baseplate, mobile tibial bearing and patella components, all of which are designed to be permanently implanted together, to achieve total reconstructive replacement of the knee joint. The BKR implant components are available in a range of sizes to accommodate anatomical variation.

The femoral and tibial baseplate components are designed for cemented fixation only, using polymethyl methacrylate (PMMA) bone cement.

Materials

Femoral and tibial baseplate components are made from cobalt chrome molybdenum alloy that complies with ISO 5832-4.

Patella and tibial bearing components are made from compression-moulded ultra high molecular weight polyethylene (UHMWPE) that complies with ISO 5834-2.

Indications

The BKR is indicated for use in patients with severe knee pain and disability due to:

- degenerative osteoarthritis, rheumatoid arthritis, or post-traumatic arthritis;
- post-traumatic loss of knee joint configuration and function; or,
- moderate valgus, varus, or flexion deformities of the knee.

Contraindications

- Previous or current history of infection in the affected joint, and/or local or systemic infection which may affect the prosthetic joint.
- Revision of failed previous total knee arthroplasty.
- Insufficient bone stock on femoral or tibial surfaces.
- Skeletal immaturity.
- Neuropathic arthropathy.
- Osteoporosis or any loss of musculature or neuromuscular disorder that compromises the affected limb.
- Severe instability of the knee joint secondary to the absence of collateral ligament integrity and function.
- Obesity - an excessively overweight or obese patient can impose forces on the prosthesis which can lead to failure of implant fixation, or failure of the implant itself.

Warnings

- The BKR implants are for single use only and must not be reused. Reuse presents a risk of

infection and a risk that implant performance may be adversely affected.

- Never reinsert an articular component that has been inserted previously, as its performance may be adversely affected.
- Polished bearing surfaces must not come into contact with hard or abrasive surfaces.
- Bearing surfaces must be clean and free of cement, or other debris prior to assembly.
- Avoid notching, scratching, or striking the implant components.
- Do not implant the femoral or tibial baseplate components without bone cement.
- Do not use the BKR implants:
 - for purposes other than the labelled indications;
 - if damage is found or caused during setup or insertion; or,
 - in conjunction with components from other knee systems.
- Avoid improper positioning or alignment of the BKR implants. The risk of implant failure is higher with inaccurate implant positioning or alignment due to unusual stress conditions which may occur, leading to a reduction in the life of the implant. Refer to the surgical technique guide for information specific to the recommended positioning and alignment of the BKR implants.

Precautions

- Before clinical use the surgeon must thoroughly understand all aspects of the surgical procedure and the limitations of the BKR implants.
- Patients should be advised that the useable life of the BKR implants is not indefinite and can be affected by various biological, mechanical and physiological factors.
- Patients should be instructed in the limitations of the implants, including potential adverse effects on the implants caused by excessive loading, as a result of patient weight or activity. Patients should be instructed to govern their activities accordingly.
- Use only surgical instruments from the dedicated BKR instrument set, as these have been specifically designed to aid in accurate insertion of the BKR implants.
- Selection of the appropriate size of implant components to use for a patient is the responsibility of the surgeon.
- Care should be taken to ensure appropriate size matching of BKR implant components. An implant size matching chart is provided to assist the user. Use of mismatched implant components may lead to poor surface contact and could result in pain, increased wear, implant instability, loss of joint function and reduced implant life.
- Selection of the thickness of the mobile tibial bearing is at the surgeons discretion - thicker bearing components may be required if the patient is young, heavy and/or physically active.

Possible Adverse Effects

- Loosening and/or disassembly of implants.
- Fracture of / damage to the implants.
- Soft tissue impingement or damage.
- Joint instability and/or dislocation.
- Malalignment of the implants.
- Bone fracture.
- Nerve damage.

- Infection.
- Soft tissue swelling or inflammation.
- Leg length discrepancies.
- Knee joint stiffness and/or limited range of motion.
- Delayed wound healing.
- Temporary or permanent neuropathy.
- Pain.
- Cardiovascular disorders including venous thrombosis, pulmonary embolism or myocardial infarction.
- Metal sensitivity.
- Corrosion of metal components.
- Osteolysis, caused by wear debris from polyethylene components.

Sterility and Shelf Life

The BKR implants are supplied in a sterile condition (sterilised by gamma irradiation) and have a shelf life of 5 years when stored unopened, under ambient conditions.

Implant packages should be inspected prior to use. Do not use an implant if the sterile packaging is open or damaged, or if the indicated expiry date has been exceeded.

Metal femoral and tibial baseplate implant components that have been opened but have not been used may be re-sterilised once (the implants must be returned to JointMedica for re-sterilisation). Patella and mobile tibial bearing implant components composed of UHMWPE must not be re-sterilised.

Manufacturer

JointMedica Limited
Unit 3 - Ball Mill Top Business Park
Main Road
Hallow
Worcestershire
WR2 6LS
United Kingdom

Tel: 01905 640008
Fax: 01905 676464
www.jointmedica.co.uk

Instructions for use cat. No.904001
Issue 11 Jun. 2019

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